

CPR/AED DRILL PROGRAM Unconscious Victim Medical Emergency

Chain of Survival

- Call 911
- Early CPR
- Early Defibrillation
- Early Advance Life Support



MEDICAL EMERGENCY DRILL AED/CPR USE

This drill is designed for any individual lay rescuer (non-medical) to practice emergency response procedures and to evaluate and improve upon the effectiveness of those procedures.

A recent Public Access Defibrillation clinical study conducted by the National Heart, Lung and Blood Institute concluded that trained laypersons can use AEDs safely and effectively. **The study resulted in twice as many people surviving out-of-hospital cardiac arrest when an AED and CPR were used over CPR alone.** Participants in this study had effective response plans that included mock cardiac arrest drills.

In 2004, the AHA's Publication – The Emergency Response Plan for Schools, emphasized routine AED practice drills and evaluations. Considering the importance and effectiveness of fire drills, a situation most anyone can relate to – the need for routine practice of rescuing a sudden cardiac arrest victim during a life or death situation is also of importance.

AED drills should be implemented by your AED coordinator and include your Building Emergency Response Team and/or Indentified First Responders. Although the main components for AED drills are standard, each building is unique and drills should be tailored to meet the needs of each facility.

MEDICAL EMERGENCY DRILL COMPONENTS:

Equipment:

- Training Manikin (preferably with clothing)
- AED Training Unit with Fast Response Kit, First Aid Kit w/ CPR face shield
- Stop watch

Communications:

- Inform main office of impending drill
- Simulate 911 call only. Do not call 911 unless prearranged with your 911 system
- Have participants state "Participating in an AED Drill" to non-participants.

Preparations:

- Training Manikins must be thoroughly cleaned for rescue breathing
- The AED coordinator should carry the AED Training Unit/Fast Response Kit

Drill Actions:

- Have rescuer retrieving the AED *place an "AED Drill in Progress" sign* in the AED storage location
- Exchange the retrieved AED with the Training Unit for rescue before using on manikin

Sample AED Drill Scenarios and Progressions:

Drills that resemble real life situations more accurately reflect the effectiveness and capabilities of the participants and the related procedures. The following are five separate drill levels. *The goal is to reach and consistently practice the Level 5 Drill.*

Drill (Level 1):

- 1. Inform selected individuals that they will be participating in an AED drill
- 2. Lead group to the drill location where you have placed a manikin
- 3. Observe the group's reactions and responses
- 4. Suggest recommendations to the rescuers
- 5. Run additional drills with multiple groups to further awareness and practice
- 6. Monitor and evaluate using Procedure Checklist and Time Line
- 7. Review, evaluate and discuss checklist / time line results and current procedures
- 8. Communicate drill results and any procedural changes with entire staff

Drill (Level 2):

Eliminate step #4 from Level 1: Monitor and evaluate using Procedure Checklist and Time Line

Drill (Level 3):

Eliminate step #1 from Level 2: Monitor and evaluate using Procedure Checklist and Time Line

Drill (Level 4):

Eliminate step #2 from Level 3 Drill: Place manikin at drill location informing first individual near scene regarding the drill. Monitor and evaluate using Procedure Checklist and Time Line

Drill (Level 5):

Same as Level 4 and include special situations such as water soaked clothing and chest; metals on manikin, simulated electrical (hot) wire cord near patient, medicine patch on chest, and/or collapse on metal bleachers. Monitor and evaluate using Procedure Checklist and Time Line

MEDICAL EMERGENCY DRILL AED/CPR USE PROCEDURE CHECK LIST AND TIME LINE		
Date Time Location		Time Line
Patient Collapses		
First Person arrives at the scene (may be first responder - not EMS):		Start Clock!
Concern for own safety considered?	Yes No	
Patient checked for responsiveness?	Yes No	
Internal call for help in accordance with emergency protocol	Yes No	
"Call 911" command given?	Yes No	
Command given to obtain AED?	Yes No	
Command given to contact First Responder?	Yes No	
Time of 911 Call		
Individual sends someone for help?	Yes No	
Individual instructs someone to meet EMS?	Yes No	
Documented emergency protocols followed?	Yes No	
Time of First Responder arrival at scene		
Concern for own safety considered?	Yes No	
Patient checked for responsiveness?	Yes No	
"Call 911" command given?	Yes No	
Time of AED command		
Patient responsiveness and breathing checked?	Yes No	
CPR started and performed correctly?	Yes No	
Documented emergency procedures followed? Time of AED arrival at Scene	Yes No	
	Vee Ne	
Clothing properly removed? Electrodes properly placed?	Yes No Yes No	
AED voice prompts followed? (especially do not touch patient)	Yes No	
Time of first AED shock	165 110	
AED voice prompts continued to be followed?	Yes No	
Patient placed in recovery position?	Yes No	
CPR performed correctly?	Yes No	
Was AED left on?	Yes No	
Monitoring of patient continued?	Yes No	
Time of EMS Arrival (Add six (6) minutes to time of 911 call)		
Were details of event properly conveyed to EMS?	Yes No	
SHOCK WITHIN 3 MINUTES / EMS WITHIN 7 MINUTES	EXCELLENT	
SHOCK WITHIN 5 MINUTES / EMS WITHIN 8 MINUTES	FAIR	
SHOCK AFTER 5 MINUTES / EMS AFTER 8 MINUTES	POOR	
YES ANSWERS 20-24		Iner Horo
YES ANSWERS 20-24 YES ANSWERS 17-19	First Voice Super Hero First Voice Sidekick	
YES ANSWERS 17-19 YES ANSWERS 12-17	First Voice Sidekick	
YES ANSWERS 6-12	You're Getting Better!	
YES ANSWERS LESS THAN 6	Keep Practici	
Other/Facility-Specific Medical Emergency Protocol Check List Addend		
	Yes No	
	Yes No	
	Yes No	
For further assistance on CPR / AED Drill Program please visit		e.us

Drill Evaluations & Comments:

Products Needed For Next Drill (reorders/restock):

Disclaimer:

This CPR / AED Drill Program is not a substitute for CPR/AED Certified Training. This is to help supplement your AED program. It is not endorsed by the American Heart Association or the American Red Cross.

For further assistance on AED Drills or AED Programs contact us at: www.firstvoice.us