



<b>Product Description</b>	Product is a medical transport device that allows for proper storage and transport of severed digits and avulsions. The device consists of a storage compartment for a cold compress (which is included with the product and sealed) and another storage compartment for gauze, saline solution, and an instruction card for proper preparation and storage of the amputated digit or avulsion for transport. Instruction card allows for recording time of injury and the injured party's name.
<b>Product Benefits</b>	<ul style="list-style-type: none"> <li>• Improves chances of reattachment by surgeon.</li> <li>• Eliminates freezing or improper storage of digits/avulsions while waiting for reattachment.</li> <li>• Decreases confusion during a multiple amputation event involving more than one victim.</li> </ul>
<b>Tamper-proof</b>	Yes. Pre-sealed package allows for a pinch-and-tear method of opening, showing package tampering or use.
<b>Dimensions</b>	Height 8 3/8" x Width 7.5" x Thickness 1"
<b>Weight</b>	11 ounces
<b>Contents</b>	Cold Compress with ammonium nitrate active ingredient (1); Gauze pads (2); Saline solution (1); Instruction card (1)
<b>Storage</b>	Preferred 0° to 50° Celsius (32° to 122° Fahrenheit)
<b>Operating Temperature</b>	Reaches 40° F and stays 50° F or under for in excess of 1 hour Refrigeration effect of ice pack lasts for up to 3 hours depending on temperature of environment surrounding use. AmpuSave™ can be put in cold ice water or ice after it has been activated to help sustain preferred operating temperatures (be careful to avoid freezing the digit/avulsion).
<b>Conditions to Avoid</b>	Heat (store in less than 50° Celsius or 122° Fahrenheit) Fire (Active ingredient is ammonium nitrate)
<b>Shelf Life</b>	3 years; Saline solution has an expiration date
<b>Testing</b>	Iowa State University validation studies can be provided upon request
<b>MSDS Information</b>	Can be provided upon request

## Accessories and Upgrades



### Celox Blood Clotting Solution

The best tool for stopping all kinds of bleeding situations.



### Emergency Instruction Device

Effectively coach users on hundreds of first aid emergencies from start to finish with the new, 1.25" thin EID.



### CPREzy CPR Assistance Device

Ensures accurate CPR pacing and depth for real-time events or training. Just turn it on, place it to the chest as suggested and begin compressing.



### Heartsine samaritan® AED Package

Be completely prepared to respond with an AED by having all the necessities in one complete package.

#### More Options

- Basic Rollout SET System
- Duffel SET System
- Rugged SET System
- Training Scenarios Book
- Waterproof EID Storage Bags
- EID Wall Mount
- EID Theft Stopper
- Custom First Aid Kits**



### Basic Bloodborne Pathogen Kit

Be prepared for the unexpected with this simple, yet powerful kit.



### First Aid Cube

With over 30 items in the Cube (more than an average kit), the First Voice® First Aid Cube is the perfect portable fit for work, sports, home, and life.

